


AGUSTA AW109/119

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.




System Overview:

- Quick change capability can be installed or removed from aircraft in approximately 15 minutes and requires minimal modifications to aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific

System Features & Options:

- Base module with pivot for easy loading
- Base is installed longitudinally
- Custom stretcher designed for patient comfort
- Infant Transport System available
- Systems Box (Base) includes Full O₂ Bottle/Equipment drawers
- LH and RH valences available



1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: Aeromed@Spectrum-Aeromed.com • www.Spectrum-Aeromed.com