

HAWKER 1000 DUAL DEDICATED

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.

System Overview:


- A dedicated configuration to specific needs and specific auxiliary medical equipment
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Allows for a variety of custom medical mounts and patient loading accessories

System Features & Options:

- Dedicated base and overhead panel
 - 7000 liters of oxygen, dual air compressors, vacuum pump and 1000 watt inverter
 - AC outlets, DC outlet and dual pneumatic outlets
 - Multiple storage drawers
- Custom stretcher specifically designed for patient comfort
- Multiple IV and medical equipment mounting options on the overhead panels
- 4 drawer storage cabinet



Dual dedicated configuration



4 drawer cabinet

1815 23rd Avenue North • Fargo, ND 58102

+1 701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: Aeromed@Spectrum-Aeromed.com • www.Spectrum-Aeromed.com