

# ITS DECK, STRETCHER BRIDGE, EQUIPMENT TABLE

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.

## Infant Transport System (ITS) Deck:

- The deck has provisions to mount an incubator and a housing for cylinders
- STC approval for various fixed wing and rotor wing aircraft
- High quality system with cost-effective installation
- Quick change capability
- The ITS deck without oxygen cylinders and incubator weighs up to 33 lbs
- Built aircraft-specific and customer specific.

## Stretcher Bridge:

- Easily attaches to Spectrum Aeromed stretcher
- Useful for mounting carry-on medical equipment
- Can hold up to 42lbs of equipment in a fixed wing aircraft
- STC approved

## Equipment Table:

- Designed with 3 belts for securing equipment
- Can securely hold up to 25lbs
- Attaches to our 2800 Module
- Weighs only 5lbs



Infant Transport System (ITS) Deck



Equipment Table (above) Stretcher Bridge (below)



1815 23rd Avenue North • Fargo, ND 58102

+1 701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: [Aeromed@Spectrum-Aeromed.com](mailto:Aeromed@Spectrum-Aeromed.com) • [www.Spectrum-Aeromed.com](http://www.Spectrum-Aeromed.com)