

STRETCHER BRIDGE

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.

System Benefits Include:

- Easily attaches to Spectrum Aeromed stretcher
- Useful for mounting carry-on medical equipment
- Can hold up to 40lbs of equipment in a fixed wing aircraft and up to 20 lbs in rotorcraft
- STC approved and data provided



Stretcher Bridge Installations



1815 23rd Avenue North • Fargo, ND 58102

+1 701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: Aeromed@Spectrum-Aeromed.com • www.Spectrum-Aeromed.com