


# AGUSTA AW139

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.



Pivot base installed longitudinally and laterally in the aircraft.

## System Overview:

- Quick change capability can be installed or removed from aircraft in approx. 15 minutes and requires minimal modifications to the aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific

## System Features & Options:

- 3 different track lengths available: 1/3 length; 1/2 length and full length
- Base module with drawers has the option of being positioned along the seat track anywhere it fits
- Base can be installed longitudinally or laterally in the aircraft
- Custom stretcher designed for patient comfort
- Infant Transport System available
- Custom equipment cabinets



1815 23rd Avenue North • Fargo, ND 58102

+1 701.235.5478 • Toll Free: 1.800.753.4340 • Fax: +1 701.235.5486

Email: [Aeromed@Spectrum-Aeromed.com](mailto:Aeromed@Spectrum-Aeromed.com) • [www.Spectrum-Aeromed.com](http://www.Spectrum-Aeromed.com)



**SPECTRUM**  
Aeromed > Supporting life.