



KING AIR 350ER

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.



KING AIR 350 Installation Configurations



System Overview:

- Quick change capability can be installed or removed from aircraft in approx. 15 minutes and requires minimal modifications to the aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Approved for up to 3 patient transport
- Approved for 2800, 20 and 2200 Series Module Systems
- Fully self-contained and lightweight with: Inverters; air and vacuum pumps; 3500 liters of oxygen and AC and DC outlets

System Features & Options:

- Med Wall available
- Incubator Transport Systems available
- 4-drawer cabinet available
- Patient loading accessories
- Stretcher bridge available



Manual loader on KA350



4-drawer cabinet installed



Dual Installation Right Side

1815 23rd Avenue North • Fargo, ND 58102

+1 701.235.5478 • Toll Free: 1.800.753.4340 • Fax: +1 701.235.5486

Email: Aeromed@Spectrum-Aeromed.com • www.Spectrum-Aeromed.com

