



# KING AIR 200

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.



KING AIR 200 Installation Configurations



## System Overview:

- Quick change capability can be installed or removed from aircraft in less than 15 minutes and requires minimal modifications to the aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Allows for a variety of custom medical mounts and patient loading accessories


## System Features & Options:

- Fully self-contained and lightweight with: 3500 liters of oxygen; dual air compressors; vacuum pump; 1000 watt inverter; AC and DC outlets and dual pneumatic outlets
- Custom stretcher specifically designed for patient comfort
- Approved for 2800 and 20 Series Module Systems
- Multiple IV Mounts

*\*Contact a Sales Representative for further Install Details*



Manual loader on KA200



Dual Installation Right Side

1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: [Aeromed@Spectrum-Aeromed.com](mailto:Aeromed@Spectrum-Aeromed.com) • [www.Spectrum-Aeromed.com](http://www.Spectrum-Aeromed.com)



**SPECTRUM**  
Aeromed > Supporting life.