

ITS DECK

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.


The Incubator Transport System (ITS) is an essential addition to any air ambulance that will be caring for infants. This incubator is a secure unit which can be installed onto a number of our modules. It allows for gentle care and ensured safety for all infant patients during those critical moments.

Infant Transport System (ITS) Deck:

- The deck has provisions to mount an incubator and cylinder housing for oxygen and air
- STC approval for various fixed wing and rotor wing
- High quality system with cost-effective installation
- Quick change capability
- The ITS deck without oxygen cylinders and incubator weighs up to 33lbs
- Built aircraft-specific and customer-specific



Infant Transport System (ITS) Deck



1815 23rd Avenue North • Fargo, ND 58102

+1 701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: Aeromed@Spectrum-Aeromed.com • www.Spectrum-Aeromed.com



SPECTRUM
Aeromed > Supporting life.