


KING AIR 200, 300

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.



KING AIR 350 Installation Configurations



System Overview:

- Quick change capability can be installed or removed from aircraft in less than 15 minutes and requires minimal modifications to the aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Allows for a variety of custom medical mounts and patient loading accessories

System Features & Options:

- Fully self-contained and lightweight
 - 3500 liters of oxygen, dual air compressors, vacuum pump and 1000 watt inverter
 - AC outlets, DC outlet and dual pneumatic outlets
- Custom stretcher specifically designed for patient comfort
- Approved for 2800 and 20 Series Module Systems
- Multiple IV Mounts



Manual Loader on KA350



Dual Installation Right Side

1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: Aeromed@Spectrum-Aeromed.com • www.Spectrum-Aeromed.com



SPECTRUM
Aeromed > Supporting life.