


# CESSNA CJ

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.



## System Overview:

- Bench Length: 48" (121.92 cm)
- Width: 17" (43.25 cm)
- Height: 10" (25.40 cm) to top of bench
- Base Unit Weight: 139 lbs (63.05 kgs)
- Outlets for electrical supply, oxygen, medical air and vacuum
- IV Mount (1)

## Spectrum Stretcher:

- Stretcher Length: 72" (183 cm.)
- Stretcher Width: 18.5" (47 cm.)
- Custom stretcher specifically designed for patient comfort with expandable armrests and pneumatically assisted fully adjustable backrest

## Electrical Supply:

- Inverter: 1000 watt
- 4 -single 230 or 115 volt outlets
- 1-28 volt DC outlet

## Seat Rail Adapter:

- Utilizes pins and fits directly into existing hardpoints on aircraft
- Aircraft specific
- Quick release capabilities 2 Year Warranty on Spectrum Aeromed parts and labor



Stretcher and base installed



1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: [Aeromed@Spectrum-Aeromed.com](mailto:Aeromed@Spectrum-Aeromed.com) • [www.Spectrum-Aeromed.com](http://www.Spectrum-Aeromed.com)



**SPECTRUM**  
Aeromed > Supporting life.