



AGUSTA AW109SP

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.



Stretcher and Base Installed



Storage options

System Overview:

- High quality, durable design and cost-effective installation
- Single and Dual Patient applications
- Built aircraft-specific and customer-specific
- Easy to install

System Features & Options:

- Folding secondary base and stretcher available; both fold for ease of storage
- Oxygen box for 2180 or 3275 Liter O2 cylinder housed in the aft cargo area
- Swiveling and Tracking Medical Attendant Seat
- 1200 Watt inverter to supply 110 volt AC Power
- 3 air compressors
- Various medical equipment mounts available



Main Base and Stretcher and Secondary Base and Stretcher Installed with Medical Attendant Seat

1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: Aeromed@Spectrum-Aeromed.com • www.Spectrum-Aeromed.com



SPECTRUM
Aeromed > Supporting life.