

AS-350 A-FRAME

A fully modular stretcher system specifically designed for the Aerospatiale AS-350 Helicopter.

Components of the A-Frame are lightweight, easy to install and fully FAA certified. With the Deluxe Spectrum Stretcher the backrest is pneumatically raised and lowered.

Once the co-pilot seat is out and the rear seat raised, all it takes is four (4) fast pins to complete the installation. The stretcher slides in and is quickly secured with an elegant fold down locking system.

Key Features:

- Deluxe Spectrum Stretcher Standard
- One Piece AS-350 Floor Adapter
- Easily Removed or Installed in less than Five Minutes
- No Structural Modification to Aircraft Required

Specifications:

- Weight: Adapter+Stretcher = 54 lb., 24.5kg
- Hard Coat Anodized Floor Adapter
- Spectrum Stretcher (17" x 72")
- Quick-Connect Hands Free Locking System
- FAA Approved Installation



1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: Aeromed@Spectrum-Aeromed.com • www.Spectrum-Aeromed.com