


SIKORSKY S76/S76C with Pivot

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.



System Overview:

- Quick change capability (20 and 2200 systems) can be installed or removed from aircraft in approximately 15 minutes and requires minimal modification to aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Can be installed with no electrical modification

System Features & Options:

- Fully self-contained and lightweight – oxygen, vacuum and air pump, AC and DC outlets and pneumatic outlets
- Custom Stretcher designed for patient comfort
- 4-drawer equipment cabinet available
- Both Med Wall Series and ITS Deck
- Base with pivot attached and additional 3-drawers for storage

Pivot Stretcher and 3 Drawer Base



1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: Aeromed@Spectrum-Aeromed.com • www.Spectrum-Aeromed.com