



# CESSNA 208, 208B

Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.



Cessna 208B



Cessna 208

## System Overview:

- Quick change capability can be installed or removed from aircraft in less than 15 minutes and requires minimal modifications to the aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific

## System Features & Options:

- Fully self-contained and lightweight
  - 3500 liters of oxygen, dual air compressors, vacuum pump and 1000 watt inverter
  - AC outlets, DC outlet and dual pneumatic outlets
- Custom stretcher specifically designed for patient comfort
- Multiple IV Mounts
- 2800, 20 and 2200 series modules available for installation
- 4 drawer equipment cabinet available
- Both Med Wall Series available
- Patient loading system available



Cessna with Short Box System (20 series) installed



2800 Series Module Installed in 208B



Patient being loaded into Caravan using Manual Patient Loader

1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: [Aeromed@Spectrum-Aeromed.com](mailto:Aeromed@Spectrum-Aeromed.com) • [www.Spectrum-Aeromed.com](http://www.Spectrum-Aeromed.com)