

INFINITY 5000X

WELCOME TO THE FUTURE OF AIR MEDICAL TECHNOLOGY

THE INFINITY 5000X

The Infinity 5000X will revolutionize the air medical industry as one of the first products designed with and built for flight nurses and air medical teams.

Since 1991, Spectrum Aeromed has engineered equipment to meet the evolving needs of air medical providers. The Infinity 5000X continues that legacy with the same dedication and high-level insight trusted worldwide.



SYSTEM OVERVIEW

- Fully self-contained with the same high-quality durable design as the previous systems.
- Quick-change configuration capability for seamless transitions between mission roles.
- Newly designed stretcher bridge allows unobstructed access when caring for patients and when transferring the patient onto and off the stretcher.
- Stretcher is specifically designed for patient comfort.
- Stretcher attaches to the base via a locking mechanism in the middle of the base system, where all controls and gauges are centrally located.
- Dual USB ports and three individual universal electrical outlets.
- 3500 liters of oxygen, dual air compressors, single vacuum pump, and 1000 watt inverter.
- Multiple IV Pole Mount locations.
- Capability to replace the stretcher with the Infant Transport System deck on the base.

With a focus on improved patient and installer experiences, Spectrum Aeromed designed the Infinity 5000X with more clearance, easier access for maintenance, and improved aesthetics.



1815 23RD AVENUE NORTH • FARGO, ND 58102 USA

+1.701.235.5478 • 1.800.753.4340

AEROMED@SPECTRUM-AEROMED.COM • SPECTRUM-AEROMED.COM